





## COMPACT WITH QUALITY PERFORMANCE

- Compact Bench Top Multi-Purpose Analyser.
- Clinical Chemistry Turbidimetry, DOA, TDM.
- Up to 360 tests/hour. (Average 300 tests/hour).
- Auto-rerun with auto-dilution.
- Stat Mode (Immediate Access- 2 min to Result).
- Unrivalled Precision CV <1%.
- Sample pre-dilution.
- Multiple sampling from each position.
- Reusable Cuvettes.
- Cost Efficient! 1.6ml Wash Solution per Test.
- Intelligent, User Friendly Software.
- Programmable Reagent Sample Mixing.
- Primary tubes 5,7,10 ml.
- Reagent Cooling - 15°C Below Ambient.  
Optional On/Off Key.
- 60 Reagent Positions (30 x 2 Reagent Sets).
- Continuous Sample Loading.
- 60 Samples Per Tray.
- Multiple Reagent Tray Programming.
- 15 Minute Incubation Time for Insensitive Assays  
i.e. ACE and Fruct.
- Vertical obstruction (clot detection).



### • HEADQUARTERS IRELAND

Audit Diagnostics,  
Business & Technology Park,  
Carrigtwohill,  
Co.Cork,  
Ireland.  
Tel +353 21 4533652  
Fax +353 21 4533653  
Email: info@auditdiagnostics.ie  
Website: www.auditdiagnostics.com

### • UK

Audit Diagnostics UK  
P.O.Box 679,  
Wivelsfield Green,  
Haywards Heath,  
West Sussex,  
RH 16 9BY,  
UK.  
Tel / Fax + 44 1444 471 494  
Email: Selwyn\_morgan@email.msn.com

### • DUBAI

Audit Diagnostics Dubai  
Intelution LLC,  
P.O.Box 73024,  
Dubai,  
United Arab Emirates.  
Tel: +971.4.3219747  
Fax: +971.4.3219689  
M: +971.50.5989219  
Email: mlateef@auditdiagnostics.ie

### • CHINA

Audit Diagnostics China  
Building 2-2512,  
No. 9 Pufang Road  
Fengtai District,  
Beijing 100078,  
China.  
Tel / Fax: 0086-10-5807 1145  
Email: minfenglin@auditdiagnostics.ie

### • INDIA

Audit Diagnostics India PVT. Ltd.  
Chanderwati Vashisht Kunj  
1st Floor  
247, Jain Nagar,  
Karala,  
Delhi-110081,  
India.  
Tel:+91(11)25474733/65164533  
Fax:+91(11)25968849  
E.Mail: vashisht@auditdiagnosticsindia.com  
web : www.auditdiagnosticsindia.com



Audit Diagnostics



## Automated Clinical Chemistry Analyser!



Immuno-Assay  
D.O.A.  
Clinical Chemistry  
Specific Proteins  
T.D.M.



Audit Diagnostics