



Introducing ClearVet™, a High Quality Liquid Based Cytology Solution for Veterinary Laboratories

- **Simple.** Fewer steps to process test sample
- **No Special Instrumentation.** Minimize upfront cost
- **Proven.** Over 250,000 successful tests to-date
- **Cost Effective.** Low cost per test
- **Increased Cellularity.** Mono layer slide for enhanced visibility

ClearVet™ Cytology

ClearVet™ has been developed over years of perfecting a hybrid process for the human clinical laboratory. With ClearVet™, labs produce a well fixed, highly cellular, and thin-layer sample free of mechanical and procedural artifacts. ClearVet™ has multiple cytology applications including: Urine, Pleural, Spinal, Synovial, FNA, and other fluids for all animals, including canine, feline, livestock, and equine.

ClearVet™ Veterinary Lab Test Kit components include: ClearVet™ Cytology Solution, Centrifuge Tube, Transport Sample Vial (preservative), Microscope Slide and Pipette



Upgrade Your Lab to ClearVet™ Today

Call +353 21 453 3652 or visit us online at www.auditdiagnostics.ie

Audit Diagnostics, Business & Technology Park
Carrigtwohill Co Cork, Ireland

