



Audit Diagnostics

Rapid test results within 5 – 10 minutes
No refrigerator needed! Room temperature is enough!

Product Name	Cat. No.	Test/Kit	Substance	Features (Storage: 2 ~ 30°C / Shelf life: 18-24 months)
Canine Heartworm Test Kit (CHW Ag)	ADA 101	10	<i>Dirofilaria immitis</i> Antigen	Detection of specific antigen of canine <i>Dirofilaria immitis</i> within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Canine Distemper Virus Test Kit (CDV Ag)	ADA 102	10	Canine Distemper Virus Antigen	Detection of specific antigen of canine distemper virus within 10 minutes Sample: ocular discharge and nasal discharge Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Canine Parvovirus Test Kit (CPV Ag)	ADA 103	10	Canine Parvovirus Antigen	Detection of specific antigen of canine parvovirus within 10 minutes Sample: feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Canine Coronavirus Test Kit (CCV Ag)	ADA 107	10	Canine Coronavirus Antigen	Detection of specific antigen of canine coronavirus within 10 minutes Sample: feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
CCV Ag / CPV Ag Test Kit (combo test)	ADA 108	10	Canine Coronavirus Antigen / Canine Parvovirus Antigen	Detection of specific antigens of canine coronavirus & parvovirus within 10 minutes Sample: feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Leishmania Test Kit (LSH Ab)	ADA 111	10	<i>Leishmania</i> Antibody	Detection of anti- <i>Leishmania</i> antibody within 10 minutes Sample: canine / feline whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Ehrlichia Test Kit (E. canis Ab)	ADA 112	10	<i>Ehrlichia canis</i> Antibody	Detection of anti- <i>Ehrlichia canis</i> antibody within 10 minutes Sample: canine / feline whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Giardia Test Kit (GIA Ag)	ADA 113	10	<i>Giardia lamblia</i> Antigen	Detection of specific antigen of <i>Giardia lamblia</i> within 10 minutes Sample: canine / feline feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Lyme Ag Test Kit (Lyme Ag)	ADA 116	10	Borrelia burgdorferi Antibody (C6 specific)	Detection of anti-Borrelia burgdorferi antibodies within 10 minutes Sample: canine whole blood, serum or plasma Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Canine Adenovirus Ag Test Kit (CAV Ag)	ADA 117	10	Canine Adenovirus Antigen	Detection of specific antigens of canine adenovirus within 10 minutes Sample: canine ocular discharge or nasal discharge Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
CAV Ag / CDV Ag Test Kit (combo test)	ADA 118	10	Canine Adenovirus Antigen / Canine Distemper virus Antigen	Detection of specific antigens of canine adenovirus and canine distemper virus within 10 minutes Sample: canine ocular discharge or nasal discharge Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Canine Pancreatic Lipase Test Kit (combo test)	ADA 121	10	Canine Pancreatic Lipase	Determination of pancreas-specific lipase levels within 10 minutes Sample: serum Contents: test cassettes, buffer bottle and disposable droppers
Relaxin Test Kit (Pregnancy Test)	ADA 122	10	Relaxin hormone (Pregnancy Test)	Detection of Relaxin hormone within 10 minutes Sample: canine and feline serum or plasma Contents: test cassettes, buffer bottle and disposable droppers

Contact:

Audit Diagnostics
Business & Technology Park
Carrigtwohill Co Cork, Ireland

Call **+353 21 453 3652**
or visit us online at:
www.auditdiagnostics.ie





Audit Diagnostics

Rapid test results within 5 – 10 minutes
No refrigerator needed! Room temperature is enough!

Product Name	Cat. No.	Test/Kit	Substance	Features (Storage: 2 ~ 30°C / Shelf life: 18-24 months)
Feline Immunodeficiency Virus Test Kit (FIV Ab)	ADA 104	10	Feline Immunodeficiency Virus Antibody	Detection of anti-feline immunodeficiency virus antibody within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Feline Leukemia Virus Test Kit (FeLV Ag)	ADA 105	10	Feline Leukemia Virus Antigen	Detection of specific antigen of feline leukemia virus within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
FeLV / FIV Ab Test Kit (combo test)	ADA 106	10	Feline Leukemia Virus Antigen / Feline Immunodeficiency Virus Antibody	Detection of feline leukemia virus antigen and anti-feline immunodeficiency virus antibody within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Feline Heartworm Test Kit (FHW Ag)	ADA 109	10	<i>Dirofilaria immitis</i> Antigen	Detection of specific antigen of feline <i>Dirofilaria immitis</i> within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Feline Parvovirus Test Kit (FPV Ag)	ADA 110	10	Feline Parvovirus Antigen	Detection of specific antigen of feline parvovirus within 10 minutes Sample: feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Leishmania Test Kit (LSH Ab)	ADA 111	10	<i>Leishmania</i> Antibody	Detection of anti- <i>Leishmania</i> antibody within 10 minutes Sample: canine / feline whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Ehrlichia Test Kit (E.canis Ab)	ADA 112	10	<i>Ehrlichia canis</i> Antibody	Detection of anti- <i>Ehrlichia canis</i> antibody within 10 minutes Sample: canine / feline whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Giardia Test Kit (GIA Ag)	ADA 113	10	<i>Giardia lamblia</i> Antigen	Detection of specific antigen of <i>Giardia lamblia</i> within 10 minutes Sample: canine / feline feces Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
FIP Test Kit (FIP Ab)	ADA 114	10	Feline Infectious Peritonitis Virus Antibody	Detection of anti-feline infectious peritonitis virus antibodies within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Feline Toxo IgG/IgM Ab Test Kit (combo test)	ADA 115	10	<i>Toxoplasma gondii</i> IgG/IgM Antibody	Detection of anti-feline <i>Toxoplasma gondii</i> IgG/IgM antibodies within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Relaxin Test Kit (Pregnancy Test)	ADA 122	10	Relaxin hormone (Pregnancy Test)	Detection of Relaxin hormone within 10 minutes Sample: canine and feline serum or plasma Contents: test cassettes, buffer bottle and disposable droppers

Contact:

Audit Diagnostics
Business & Technology Park
Carrigtwohill Co Cork, Ireland

Call **+353 21 453 3652**
or visit us online at:
www.auditdiagnostics.ie





Audit Diagnostics

High sensitivity and specificity
Supply of long expiry dated products. More than 1 year!

Product Name	Cat. No.	Test/Kit	Substance	Features (Storage: 2 ~ 30°C / Shelf life: 18-24 months)
Avian Influenza Ag Test Kit (AI Ag)	ADLS 01	10	Avian Influenza Virus Antigen (group A)	Detection of specific antigens of avian influenza virus within 10 minutes Sample: feces or cloaca swab Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
NDV Ag Test Kit	ADLS 02	10	Newcastle disease Virus Antigen	Detection of specific antigens of newcastle disease virus within 10 minutes Sample: avian cloaca for trachea swab Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
IBDV AG Test Kit	ADLS 03	10	Infectious bursal disease virus Antigen	Detection of specific antigens of infectious bursal disease virus within 10 minutes Sample: avian bursal of fabricious or cloaca swab Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
FMDV Ag Test Kit	ADLS 04	10	Foot-and-Mouth Disease NSP antibody	Detection of anti-FMDV antibodies within 10 minutes Sample: bovine, ovine and porcine whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Swine Influenza Ag Test Kit (SI Ag)	ADLS 05	10	Swine Influenza Virus Antigen	Detection of specific antigen of Swine influenza virus within 10 minutes Sample: ocular discharge and nasal discharge Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Swine Influenza Virus type-specific Test Kit (SI Ag)	ADLS 06	10	Swine Influenza Virus Antigen (HA1, HA3, NA1, NA2)	Detection of type-specific (HA1, HA3, NA1, NA2) antigen of Swine influenza virus within 10 minutes Sample: ocular discharge and nasal discharge Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Classical Swine Fever Virus Ag Test Kit (CSFV Ab)	ADLS 07	10	Classical Swine Fever Virus Antigen	Detection of specific antigens of classical swine fever virus within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Classical Swine Fever Virus Ab Test Kit (CSFV Ab)	ADLS 08	10	Classical Swine Fever Virus Antibody	Detection of specific antibodies of classical swine fever virus within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Classical Swine Fever Virus Ab_DF Test Kit (CSFV Ab_DF)	ADLS 09	10	Classical Swine Fever Virus Antigen (titer ≥ 1:32)	Detection of specific antibodies (titer ≥ 1:32) of classical swine fever virus within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Swine Toxo Ab Test Kit	ADLS 10	10	<i>Toxoplasma gondii</i> Antibody (IgG)	Detection of anti-swine <i>Toxoplasma gondii</i> antibodies within 10 minutes Sample: whole blood, serum or plasma Contents: test cassettes, buffer bottle and disposable droppers
Avian Influenza H5 Ag Test Kit (AI H5 Ag)	ADLS 11	10	Avian Influenza H5 Virus Antigen	Detection of H5 specific antigens avian influenza virus within 10 minutes Sample: feces or cloaca swab Contents: test cassettes, buffer tubes, disposable droppers and cotton swabs
Rabies Ag Test Kit	ADA 119	10	Rabies virus Antigen	Detection of specific antigens of rabies virus within 10 minutes Sample: canine, bovine, porcine's feces Contents: test cassettes,buffer tubes, disposable droppers and cotton swabs
Rotavirus Ag Test Kit	ADA 120	10	Rotavirus Antigen	Detection of specific antigens of rotavirus within 10 minutes Sample: canine, bovine, porcine's feces Contents: test cassettes,buffer tubes, disposable droppers and cotton swabs

Contact:

Audit Diagnostics
Business & Technology Park
Carrigtwohill Co Cork, Ireland

Call **+353 21 453 3652**
or visit us online at:
www.auditdiagnostics.ie

